

Anti-CD162/Selplg antibody [4RA10] {PE/Cyanine5.5} (STJA0001410)

STJA0001410

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-P-selectin glycoprotein ligand 1 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	4RA10
Concentration	0.2 mg/mL
Conjugation	PE/Cyanine5.5
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	NA
Gene Symbol	Selplg
Uniprot ID	SELPL_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	