

Anti-CD120b/Tnfrsf1b antibody [TR75-54.7] (STJA0001413)

STJA0001413

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Armenian Hamster monoclonal antibody anti-Tumor necrosis factor receptor superfamily member 1B is suitable for use in Flow Cytometry, Western Blot, Universal and Immunoprecipitation research applications.
Applications	FC, WB, Block, IP
Host/Source	Armenian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	TR75-54.7
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	21938
Gene Symbol	Tnfrsf1b
Uniprot ID	TNR1B_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	