

## Anti-CD120b/Tnfrsf1b antibody [TR75-54.7] (STJA0001414)

STJA0001414

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Armenian Hamster monoclonal antibody anti-Tumor necrosis factor receptor superfamily member 1B is suitable for use in Flow Cytometry, Western Blot and Immunoprecipitation research applications.
<b>Applications</b>	FC, WB, IP
<b>Host/Source</b>	Armenian Hamster
<b>Reactivity</b>	Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	TR75-54.7
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution</b>	
<b>Range</b>	
<b>Formulation</b>	PBS with 0.05% Proclin300
<b>Isotype</b>	IgG
<b>Storage</b>	Recommend storing at 2-8°C. Do not freeze this product.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">21938</a>
<b>Gene Symbol</b>	<a href="#">Tnfrsf1b</a>
<b>Uniprot ID</b>	<a href="#">TNR1B_MOUSE</a>
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	