

**Anti-CD120b/Tnfrsf1b antibody [TR75-54.7] {Biotin} (STJA0001415)**

STJA0001415

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Armenian Hamster monoclonal antibody anti-Tumor necrosis factor receptor superfamily member 1B is suitable for use in Flow Cytometry research applications.
<b>Applications</b>	FC
<b>Host/Source</b>	Armenian Hamster
<b>Reactivity</b>	Mouse

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	TR75-54.7
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Biotin
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with 0.05% Proclin300, 1% BSA
<b>Isotype</b>	IgG
<b>Storage</b>	Recommend storing at 2-8°C. Do not freeze this product.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	21938
<b>Gene Symbol</b>	Tnfrsf1b
<b>Uniprot ID</b>	TNR1B_MOUSE
<b>Immunogen</b>	
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	