

Anti-CD29/ITGB1 antibody [TS2/16.2.1] {PE/Cyanine5.5} (STJA0001438)

STJA0001438

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Integrin beta-1 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	TS2/16.2.1
Concentration	
Conjugation	PE/Cyanine5.5
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	3688
Gene Symbol	ITGB1
Uniprot ID	ITB1_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	