

Anti-CD62L/SELL antibody [DREG56] {AF647} (STJA0001445)

STJA0001445

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-L-selectin is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	DREG56
Concentration	
Conjugation	AF647
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	6402
Gene Symbol	SELL
Uniprot ID	LYAM1_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	