

Anti-CD38 antibody [HIT2] {Biotin} (STJA0001467)

STJA0001467

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 1 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HIT2
Concentration	0.5 mg/mL
Conjugation	Biotin
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	952
Gene Symbol	CD38
Uniprot ID	CD38_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	