

Anti-CD47 antibody [CC2C6D4] {Biotin} (STJA0001475)

STJA0001475

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Leukocyte surface antigen CD47 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	CC2C6D4
Concentration	0.5 mg/mL
Conjugation	Biotin
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	961
Gene Symbol	CD47
Uniprot ID	CD47_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	