

Anti-CD35/CR1 antibody [E11] (STJA0001479)

STJA0001479

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Complement receptor type 1 is suitable for use in Flow Cytometry, Immunohistochemistry and Immunoprecipitation research applications.
Applications	FC, IHC, IP
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	E11
Concentration	0,5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG1k
Storage	Recommend storing at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	1378
Gene Symbol	CR1
Uniprot ID	CR1_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	