

Anti-CD35/CR1 antibody [E11] {AF647} (STJA0001484)

STJA0001484

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Complement receptor type 1 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	E11
Concentration	
Conjugation	AF647
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	1378
Gene Symbol	CR1
Uniprot ID	CR1_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	