

Anti-CD57/B3GAT1 antibody [HNK-1] (STJA0001485)

STJA0001485

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Galactosylgalactosylxylosylprotein 3-beta-glucuronosyltransferase 1 is suitable for use in Flow Cytometry, Western Blot and Immunohistochemistry research applications.
Applications	FC, WB, IHC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HNK-1
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution	
Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgMk
Storage	Recommend storing at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	27087
Gene Symbol	B3GAT1
Uniprot ID	B3GA1_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	