

Anti-CD58 antibody [TS2/9.1] (STJA0001489)

STJA0001489

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Lymphocyte function-associated antigen 3 is suitable for use in Flow Cytometry and Immunoprecipitation research applications.
Applications	FC, IP
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	TS2/9.1
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG1k
Storage	Recommend storing at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	965
Gene Symbol	CD58
Uniprot ID	LFA3_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	