

Anti-CD58 antibody [TS2/9.1] {Biotin} (STJA0001490)

STJA0001490

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Lymphocyte function-associated antigen 3 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	TS2/9.1
Concentration	0.5 mg/mL
Conjugation	Biotin
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	965
Gene Symbol	CD58
Uniprot ID	LFA3_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	