

Anti-CD74 antibody [LN2] (STJA0001493)

STJA0001493

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-HLA class II histocompatibility antigen gamma chain is suitable for use in Flow Cytometry, Immunohistochemistry, Immunofluorescence and Immunoprecipitation research applications.
Applications	FC, IHC, IF, IP
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	LN2
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG1k
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	972
Gene Symbol	CD74
Uniprot ID	HG2A_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	