

Anti-CD74 antibody [LN2] {PE/Cyanine7} (STJA0001496)

STJA0001496

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-HLA class II histocompatibility antigen gamma chain is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	LN2
Concentration	
Conjugation	PE/Cyanine7
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	972
Gene Symbol	CD74
Uniprot ID	HG2A_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	