

Anti-CD32/FCGR2B antibody [IV-3] {FITC} (STJA0001502)

STJA0001502

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Low affinity immunoglobulin gamma Fc region receptor II-b is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	IV-3
Concentration	
Conjugation	FITC
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2bk
Storage	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	2213
Gene Symbol	FCGR2B
Uniprot ID	FCG2B_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	