

Anti-CD11b antibody [M1/70] (STJA0001524)

STJA0001524

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-CD11b is suitable for use in Flow Cytometry, Immunohistochemistry and Immunoprecipitation research applications.
Applications	FC, IHC, IP
Host/Source	Rat
Reactivity	Mouse, Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M1/70
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG2bk
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	NA
Gene Symbol	Itgam
Uniprot ID	ITAM_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	