

Anti-CD41/ITGA2B antibody [HIP8] {PE/Cyanine5.5} (STJA0001540)

STJA0001540

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Integrin alpha-IIb is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HIP8
Concentration	
Conjugation	PE/Cyanine5.5
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	3674
Gene Symbol	ITGA2B
Uniprot ID	ITA2B_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	