

Anti-CD103/Itgae antibody [M290] (STJA0001551)

STJA0001551

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Integrin alpha-E is suitable for use in Flow Cytometry, Immunohistochemistry and Immunoprecipitation research applications.
Applications	FC, IHC, IP
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M290
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG2ak
Storage	Recommend storing at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	16407
Gene Symbol	Itgae
Uniprot ID	ITAE_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	