

Anti-CD106/Vcam1 antibody [M/K-2.7] (STJA0001564)

STJA0001564

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Vascular cell adhesion protein 1 is suitable for use in Flow Cytometry and Immunofluorescence research applications.
Applications	FC, IF
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M/K-2.7
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG1k
Storage	Recommend storing at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	NA
Gene Symbol	Vcam1
Uniprot ID	VCAM1_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	