

## Anti-CD117/Kit antibody [2B8] (STJA0001571)

STJA0001571

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rat monoclonal antibody anti-Mast/stem cell growth factor receptor Kit is suitable for use in Flow Cytometry, Immunohistochemistry, Immunofluorescence and Immunoprecipitation research applications.
<b>Applications</b>	FC, IHC, IF, IP
<b>Host/Source</b>	Rat
<b>Reactivity</b>	Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	2B8
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution</b>	
<b>Range</b>	
<b>Formulation</b>	PBS with 0.05% Proclin300
<b>Isotype</b>	IgG2bk
<b>Storage</b>	Recommend storing at 2-8°C. Do not freeze this product.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">16590</a>
<b>Gene Symbol</b>	<a href="#">Kit</a>
<b>Uniprot ID</b>	<a href="#">KIT_MOUSE</a>
<b>Immunogen</b>	
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	