

**Anti-CD117/Kit antibody [2B8] {AF647} (STJA0001580)**

STJA0001580

**GENERAL INFORMATION**

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Primary antibodies  |
| <b>Short Description</b> | Rat monoclonal antibody anti-Mast/stem cell growth factor receptor Kit is suitable for use in Flow Cytometry research applications. |
| <b>Applications</b>      | FC  |
| <b>Host/Source</b>       | Rat   |
| <b>Reactivity</b>        | Mouse   |

**PRODUCT PROPERTIES**

|                            |  |
|----------------------------|--|
| <b>Clonality</b>           | Monoclonal   |
| <b>Clone ID</b>            | 2B8  |
| <b>Concentration</b>       | 0.5 mg/mL  |
| <b>Conjugation</b>         | AF647  |
| <b>Purification</b>        |  |
| <b>Dilution Range</b>      |  |
| <b>Formulation</b>         | PBS with 0.05% Proclin300, 1% BSA  |
| <b>Isotype</b>             | IgG2bk   |
| <b>Storage Instruction</b> | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

**TARGET INFORMATION**

|                           |                           |
|---------------------------|---------------------------|
| <b>Gene ID</b>            | <a href="#">16590</a>     |
| <b>Gene Symbol</b>        | <a href="#">Kit</a>       |
| <b>Uniprot ID</b>         | <a href="#">KIT_MOUSE</a> |
| <b>Immunogen</b>          |                           |
| <b>Immunogen Region</b>   |                           |
| <b>Specificity</b>        |                           |
| <b>Immunogen Sequence</b> |                           |