

## Anti-CD71/Tfrc antibody [R17 217.1.3/TIB-219] (STJA0001582)

STJA0001582

### GENERAL INFORMATION

|                          |  |
|--------------------------|--|
| <b>Product Type</b>      | Primary antibodies   |
| <b>Short Description</b> | Rat monoclonal antibody anti-Transferrin receptor protein 1 is suitable for use in Flow Cytometry research applications. |
| <b>Applications</b>      | FC   |
| <b>Host/Source</b>       | Rat  |
| <b>Reactivity</b>        | Mouse  |

### PRODUCT PROPERTIES

|                            |   |
|----------------------------|---|
| <b>Clonality</b>           | Monoclonal  |
| <b>Clone ID</b>            | R17 217.1.3/TIB-219                                     |
| <b>Concentration</b>       | 0.5 mg/mL   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Purification</b>        |   |
| <b>Dilution Range</b>      |   |
| <b>Formulation</b>         | PBS with 0.05% Proclin300                               |
| <b>Isotype</b>             | IgG2ak  |
| <b>Storage Instruction</b> | Recommend storing at 2-8°C. Do not freeze this product. |

### TARGET INFORMATION

|                           |                            |
|---------------------------|----------------------------|
| <b>Gene ID</b>            | <a href="#">22042</a>      |
| <b>Gene Symbol</b>        | <a href="#">Tfrc</a>       |
| <b>Uniprot ID</b>         | <a href="#">TFR1_MOUSE</a> |
| <b>Immunogen</b>          |                            |
| <b>Immunogen Region</b>   |                            |
| <b>Specificity</b>        |                            |
| <b>Immunogen Sequence</b> |                            |