

**Anti-CD71/Tfrc antibody [R17 217.1.3/TIB-219] {Biotin} (STJA0001583)**

STJA0001583

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rat monoclonal antibody anti-Transferrin receptor protein 1 is suitable for use in Flow Cytometry research applications.
<b>Applications</b>	FC
<b>Host/Source</b>	Rat
<b>Reactivity</b>	Mouse

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	R17 217.1.3/TIB-219
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Biotin
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with 0.05% Proclin300, 1% BSA
<b>Isotype</b>	IgG2ak
<b>Storage Instruction</b>	Recommend storing at 2-8°C. Do not freeze this product.

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">22042</a>
<b>Gene Symbol</b>	<a href="#">Tfrc</a>
<b>Uniprot ID</b>	<a href="#">TFR1_MOUSE</a>
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	