

Anti-CD90.2 antibody [30H12] (STJA0001590)

STJA0001590

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-CD90.2 is suitable for use in Universal, Flow Cytometry and Immunohistochemistry research applications.
Applications	Depletion, FC, IHC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	30H12
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG2bk
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	21838
Gene Symbol	Thy1
Uniprot ID	THY1_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	