

Anti-Cd4 antibody [GK1.5] (STJA0001616)

STJA0001616

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-T-cell surface glycoprotein CD4 is suitable for use in Flow Cytometry, Western Blot, Immunohistochemistry and Immunoprecipitation research applications.
Applications	FC, WB, IHC, IP
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	GK1.5
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG2bk
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	12504
Gene Symbol	Cd4
Uniprot ID	CD4_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	