

Anti-Ifng antibody [X mg1.2] (STJA0001645)

STJA0001645

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Rat monoclonal antibody anti-Interferon gamma is suitable for use in IConfocal MicroscopyC, ELISA, Western Blot and Immunohistochemistry research applications. |
| Applications | ICFC, ELISA, WB, IHC |
| Host/Source | Rat |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Monoclonal |
| Clone ID | X mg1.2 |
| Concentration | 0.5 mg/mL |
| Conjugation | Unconjugated |
| Purification | |
| Dilution Range | |
| Formulation | PBS with 0.05% Proclin300 |
| Isotype | IgG1k |
| Storage Instruction | Recommend storing at 2-8°C. Do not freeze this product. |

TARGET INFORMATION

| | |
|---------------------------|----------------------------|
| Gene ID | 15978 |
| Gene Symbol | Ifng |
| Uniprot ID | IFNG_MOUSE |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |