

Anti-Ifng antibody [X mg1.2] {Biotin} (STJA0001646)

STJA0001646

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Interferon gamma is suitable for use in IConfocal MicroscopyC research applications.
Applications	ICFC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	X mg1.2
Concentration	0.5 mg/mL
Conjugation	Biotin
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	15978
Gene Symbol	Ifng
Uniprot ID	IFNG_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	