

Anti-CD25/Il2ra antibody [PC-61.5.3] (STJA0001653)

STJA0001653

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Interleukin-2 receptor subunit alpha is suitable for use in Universal, Universal, Flow Cytometry, Immunohistochemistry and Immunoprecipitation research applications.
Applications	Block, Depletion, FC, IHC, IP
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	PC-61.5.3
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG1k
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	16184
Gene Symbol	Il2ra
Uniprot ID	IL2RA_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	