

Anti-CD3 Epsilon antibody [145-2C11] {AF647} (STJA0001665)

STJA0001665

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Armenian Hamster monoclonal antibody anti-T-cell surface glycoprotein CD3 epsilon chain is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Armenian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	145-2C11
Concentration	
Conjugation	AF647
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12501
Gene Symbol	Cd3e
Uniprot ID	CD3E_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	