

Anti-Cd8a antibody [53-6.7] (STJA0001671)

STJA0001671

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-T-cell surface glycoprotein CD8 alpha chain is suitable for use in Flow Cytometry, Western Blot, Immunofluorescence and Immunoprecipitation research applications.
Applications	FC, WB, IF, IP
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	53-6.7
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG2ak
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	12525
Gene Symbol	Cd8a
Uniprot ID	CD8A_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	