

Anti-Cd4 antibody [OX-38] (STJA0001675)

STJA0001675

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-T-cell surface glycoprotein CD4 is suitable for use in Flow Cytometry, Immunohistochemistry and Immunoprecipitation research applications.
Applications	FC, IHC, IP
Host/Source	Mouse
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	OX-38
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG2ak
Storage	Recommend storing at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	24932
Gene Symbol	Cd4
Uniprot ID	CD4_RAT
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	