

Anti-CD115/Csf1r antibody [AFS98] (STJA0001690)

STJA0001690

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Macrophage colony-stimulating factor 1 receptor is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	AFS98
Concentration	0,5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG2ak
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	12978
Gene Symbol	Csf1r
Uniprot ID	CSF1R_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	