

Anti-Ly6g antibody [1A8] (STJA0001696)

STJA0001696

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Lymphocyte antigen 6G is suitable for use in Flow Cytometry, Immunohistochemistry and Immunofluorescence research applications.
Applications	FC, IHC, IF
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	1A8
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG2ak
Storage	Recommend storing at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	546644
Gene Symbol	Ly6g
Uniprot ID	LY6G_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	