

## Anti-HLA-DRA antibody [L243] (STJA0001708)

STJA0001708

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-HLA class II histocompatibility antigen, DR alpha chain is suitable for use in Flow Cytometry, Western Blot, Immunohistochemistry and Immunoprecipitation research applications.
<b>Applications</b>	FC, WB, IHC, IP
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	L243
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution</b>	
<b>Range</b>	
<b>Formulation</b>	PBS with 0.05% Proclin300
<b>Isotype</b>	IgG2ak
<b>Storage</b>	Recommend storing at 2-8°C. Do not freeze this product.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">3122</a>
<b>Gene Symbol</b>	<a href="#">HLA-DRA</a>
<b>Uniprot ID</b>	<a href="#">DRA_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	