

Anti-CD18/Itgb2 antibody [M18/2] (STJA0001718)

STJA0001718

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Integrin beta-2 is suitable for use in Flow Cytometry, Western Blot, Immunohistochemistry and Immunoprecipitation research applications.
Applications	FC, WB, IHC, IP
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M18/2
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG2ak
Storage	Recommend storing at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	16414
Gene Symbol	Itgb2
Uniprot ID	ITB2_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	