

Anti-CD18/Itgb2 antibody [M18/2] {Biotin} (STJA0001720)

STJA0001720

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Integrin beta-2 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M18/2
Concentration	0.5 mg/mL
Conjugation	Biotin
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2ak
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	16414
Gene Symbol	Itgb2
Uniprot ID	ITB2_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	