

## Anti-Ly-6G/Ly-6C (Gr-1) antibody [RB6-8C5] (STJA0001753)

STJA0001753

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rat monoclonal antibody anti-Ly-6G/Ly-6C (Gr-1) is suitable for use in Universal, Flow Cytometry, Western Blot, Immunohistochemistry and Immunoprecipitation research applications.
<b>Applications</b>	Depletion, FC, WB, IHC, IP
<b>Host/Source</b>	Rat
<b>Reactivity</b>	Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	RB6-8C5
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
<b>Isotype</b>	IgG2bk
<b>Storage</b>	Recommend storing at 2-8°C. Do not freeze this product.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	546644
<b>Gene Symbol</b>	Ly6g
<b>Uniprot ID</b>	LY6G_MOUSE
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	