

Anti-Ly-6G/Ly-6C (Gr-1) antibody [RB6-8C5] (STJA0001754)

STJA0001754

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rat monoclonal antibody anti-Ly-6G/Ly-6C (Gr-1) is suitable for use in Flow Cytometry, Western Blot, Immunohistochemistry and Immunoprecipitation research applications. |
| Applications | FC, WB, IHC, IP |
| Host/Source | Rat |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|-----------------------|---|
| Clonality | Monoclonal |
| Clone ID | RB6-8C5 |
| Concentration | 0.5 mg/mL |
| Conjugation | Unconjugated |
| Purification | |
| Dilution Range | |
| Formulation | PBS with 0.05% Proclin300 |
| Isotype | IgG2bk |
| Storage | Recommend storing at 2-8°C. Do not freeze this product. |
| Instruction | |

TARGET INFORMATION

| | |
|--------------------|----------------------------|
| Gene ID | 546644 |
| Gene Symbol | Ly6g |
| Uniprot ID | LY6G_MOUSE |
| Immunogen | |
| Immunogen | |
| Region | |
| Specificity | |
| Immunogen | |
| Sequence | |