

Anti-Ly-6G/Ly-6C (Gr-1) antibody [RB6-8C5] {PE/Cyanine5} (STJA0001756)

STJA0001756

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Ly-6G/Ly-6C (Gr-1) is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	RB6-8C5
Concentration	
Conjugation	PE/Cyanine5
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2bk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	546644
Gene Symbol	Ly6g
Uniprot ID	LY6G_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	