

Anti-Ly6C antibody [Monts 1] (STJA0001759)

STJA0001759

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Ly6C is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Monts 1
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG2ak
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	17067
Gene Symbol	Ly6c1
Uniprot ID	LY6C1_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	