

Anti-TCR Gamma/Delta antibody [UC7-13D5] (STJA0001778)

STJA0001778

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Armenian Hamster monoclonal antibody anti-TCR Gamma/Delta is suitable for use in Flow Cytometry and Immunoprecipitation research applications.
Applications	FC, IP
Host/Source	Armenian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	UC7-13D5
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	