

Anti-CD19 antibody [4G7] (STJA0001801)

STJA0001801

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-B-lymphocyte antigen CD19 is suitable for use in Flow Cytometry and Immunofluorescence research applications.
Applications	FC, IF
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	4G7
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG1k
Storage	Recommend storing at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	930
Gene Symbol	CD19
Uniprot ID	CD19_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	