

Anti-CD19-antibody [4G7] {PE/TR} (STJA0001805)

STJA0001805

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD19 is suitable for use in Flow Cytometry.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	4G7
Concentration	
Conjugation	PE/TR
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1, k
Storage Instruction	Store at 2-8°C and protected from prolonged exposure to light. Do not freeze.

TARGET INFORMATION

Gene ID	930
Gene Symbol	CD19
Uniprot ID	CD19_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	