

Anti-CD154/CD40LG antibody [5C8] (STJA0001806)

STJA0001806

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD40 ligand is suitable for use in Universal, Flow Cytometry and Immunoprecipitation research applications.
Applications	Block, FC, IP
Host/Source	Mouse
Reactivity	Human, Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	5C8
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG2ak
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	959
Gene Symbol	CD40LG
Uniprot ID	CD40L_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	