

Anti-HLA-A, B, C antibody [W6/32] (STJA0001812)

STJA0001812

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-HLA-A, B, C is suitable for use in Flow Cytometry, Western Blot, Immunohistochemistry and Immunoprecipitation research applications.
Applications	FC, WB, IHC, IP
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	W6/32
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG2ak
Storage	Recommend storing at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	3106
Gene Symbol	HLA-B
Uniprot ID	HLAB_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	