

Anti-CD279/Pdcd1 antibody [29F.1A12] (STJA0001817)

STJA0001817

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Programmed cell death protein 1 is suitable for use in Universal, Universal, Flow Cytometry, Western Blot, Immunohistochemistry and Immunofluorescence research applications.
Applications	Block, Neut, FC, WB, IHC, IF
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	29F.1A12
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG2ak
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	18566
Gene Symbol	Pdcd1
Uniprot ID	PDCD1_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	