

## Anti-CD279/Pdcd1 antibody [29F.1A12] (STJA0001818)

STJA0001818

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rat monoclonal antibody anti-Programmed cell death protein 1 is suitable for use in Flow Cytometry, Western Blot, Immunohistochemistry and Immunofluorescence research applications.
<b>Applications</b>	FC, WB, IHC, IF
<b>Host/Source</b>	Rat
<b>Reactivity</b>	Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	29F.1A12
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with 0.05% Proclin300
<b>Isotype</b>	IgG2ak
<b>Storage</b>	Recommend storing at 2-8°C. Do not freeze this product.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">18566</a>
<b>Gene Symbol</b>	<a href="#">Pdcd1</a>
<b>Uniprot ID</b>	<a href="#">PDCD1_MOUSE</a>
<b>Immunogen</b>	
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	