

Anti-CD279/Pdcd1 antibody [29F.1A12] {Biotin} (STJA0001819)

STJA0001819

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Programmed cell death protein 1 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	29F.1A12
Concentration	0.5 mg/mL
Conjugation	Biotin
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2ak
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	18566
Gene Symbol	Pdcd1
Uniprot ID	PDCD1_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	