

Anti-Cd274/PD-L1 antibody [10F.9G2] (STJA0001821)

STJA0001821

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Programmed cell death 1 ligand 1 is suitable for use in Flow Cytometry, Immunohistochemistry and Immunofluorescence research applications.
Applications	FC, IHC, IF
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	10F.9G2
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG2bk
Storage	Recommend storing at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	60533
Gene Symbol	Cd274
Uniprot ID	PD1L1_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	